


COMP	COMPONENT NAME: MASTERFLEX® L/S® Y-TUBING ASSEMBLY								
Item Key	Part Description	Part #	Total Qty in Assem bly	Length of Tubing	Notes:				
1	Tubing, Masterflex L/S		2	17.0" (+1.0"/ -0.0")					
2	Tubing, Masterflex L/S		2	32.0" (+1.0"/ -0.0")					
3	Y-Fitting		2		Y Fitting with appropriate banding to ensure engagement of tubing to the fitting.				
	Overall the formulations that we would be interested in would be: PharMed BPT and BioPharm Plus								

Assembly Notes: 1. Attach tubing segments to Y-fitting inserts as shown. 2. Complete with appropriate banding to ensure a leakproof connection when under tubing rated pressures.



COMPONENT NAME: MASTERFLEX® I/P® Y-TUBING ASSEMBLY								
Item Key	Part Description	Part #	Total Qty in Assem bly	Length of Tubing	Notes:			
1	Tubing, Masterflex I/P		2	17.0" (+1.0"/ -0.0")				
2	Tubing, Masterflex I/P		2	32.0" (+1.0"/ -0.0")				
3	Y - Fitting		2		Y Fitting with appropriate banding to ensure engagement of tubing to the fitting.			
	Overall the formulations that we would be interested in would be: PharMed BPT and BioPharm Plus							
	PharMed BPT and							

Assembly Notes: Assembly Notes: 1. Attach tubing segments to Y-fitting inserts as shown. 2. Complete with appropriate banding to ensure a leakproof connection when under tubing rated pressures.